

Technical Data

Product: GRAPHIX ONE
Part No.: 230680V01

Display range:	10 ⁻¹⁰ - 2000 mbar
Number of measurement channels:	1
Possible active sensors:	CTR 90/91, CTR 100/101 series DU 200/201, DU2000/2001/2001 rel. Series TTR 90/91/96, TTR 100/101 series TTR 211/216/, TTR 911/916 series PTR 90, PTR 225 series ITR 90/200 series
Unit of measurement:	mbar, Torr, Micron, Pascal, Psi
Interfaces:	RS 232, RS 485
Relay interface:	25 pol. Sub-D Buchse
Relay load rating:	1A, 30V AC/30V DC
Mains connection:	100-240V, 50/60Hz
Protection class:	IP 20
Dimensions: (W x H x D)	106,4 x 128,5 x 164,5 mm

Technical data are subject to change